
## Polycheck® Allergy Diagnostics Multiparameter Technology

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Inhalation 30-I 

t23	Cypress pollen
t04	Hazel pollen
t15	White Ash pollen
t07	White Oak pollen
t09	Olive pollen
t03	Birch pollen
g12	Rye pollen
g14	Cultivated Oat
g08	June Grass pollen
g06	Timothy Grass pollen
g15	Cultivated Wheat pollen
g03	Cocksfoot pollen
g02	Bermuda Grass pollen
w09	Plantain pollen
w10	Goosefoot pollen
w21	Wall Pellitory pollen
w01	Ragweed pollen
w06	Mugwort pollen
e01	Cat epithelia
e02/e05	Dog epithelia
e03	Horse epithelia
m06	<i>Alternaria alternata</i>
m03	<i>Aspergillus fumigatus</i>
m05	<i>Candida albicans</i>
d201	<i>Blomia tropicalis</i>
d01	<i>D. pteronyssinus</i>
d02	<i>D. farinae</i>
k82	Latex
i06	Cockroach
k202	CCD <sup>9</sup>

Inhalation 10-I 

t03	Birch pollen
t02	Alder pollen
t04	Hazel pollen
t07	White Oak pollen
t09	Olive pollen
g06	Timothy Grass pollen
g12	Rye pollen
w01	Ragweed pollen
w06	Mugwort pollen
w09	Plantain pollen

Inhalation 10-II 

d01	<i>D. pteronyssinus</i>
d02	<i>D. farinae</i>
e02/e05	Dog epithelia
e01	Cat epithelia
e03	Horse epithelia
e81	Sheep epithelia
m03	<i>Aspergillus fumigatus</i>
m02	<i>Cladosporium herbarum</i>
m01	<i>Penicillium notatum</i>
m06	<i>Alternaria alternata</i>

Inhalation 10-III 

t03	Birch pollen
gx7	6 Grass Mix <sup>1</sup>
g12	Rye pollen
w06	Mugwort pollen
d01	<i>D. pteronyssinus</i>
d02	<i>D. farinae</i>
e02/e05	Dog epithelia
e01	Cat epithelia
ex78	Feathers Mix <sup>2</sup>
m02/m06	<i>Cladosporium herbarum</i> / <i>Alternaria alternata</i> Mix

Food 20-I 

f17	Hazelnut
f13	Peanut
f16	Walnut
f20	Almond
f02	Milk
f01	Egg White
f75	Egg Yolk
f78	Casein
f35	Potato
f85	Celery
f31	Carrot
f25	Tomato
f03	Codfish
f24	Shrimp
f95	Peach
f49	Apple
f14	Soybean
f04	Wheat flour
f10	Sesame
f05	Rye flour

Food 10-III 

f02	Milk
f01	Egg White
f75	Egg Yolk
f78	Casein
f14	Soya
f09	Rice
f93	Cacao [Chocolate]
f49	Apple
f31	Carrot
fx10	Flour Mix <sup>4</sup>

Food 10-IV 

f13	Peanut
f10	Sesame
f26	Pork meat
f83	Chicken meat
f27	Beef meat
f03	Codfish
fx29	Citrus Mix <sup>3</sup>
f85	Celery
f95	Peach
f25	Tomato

## Food 30-I

f25	Tomato
f96	Avocado
f92	Banana
fx29	Citrus Mix <sup>3</sup>
f84	Kiwi
f13	Peanut
f17	Hazelnut
f12	Pea
f14	Soybean
f85	Celery
f27	Beef meat
f83	Chicken meat
f26	Pork meat
f207	Clam Shell
f24	Shrimp
f23	Crab
f40	Tuna, Yellowfin
f03	Codfish
f47	Garlic
f48	Onion
f45	Baker's Yeast
f10	Sesame
f09	Rice
f08	Corn
f04	Wheat flour
f76	Alpha-Lactalbumin
f77	Beta-Lactoglobulin
f78	Casein
f01	Egg White
k202	CCD <sup>9</sup>

## Paediatric 20-I

f13	Peanut
f02	Milk
f01	Egg White
f75	Egg Yolk
f35	Potato
f31	Carrot
f03	Codfish
f49	Apple
f14	Soybean
f04	Wheat flour
t03	Birch pollen
g06	Timothy Grass pollen
w06	Mugwort pollen
d01	D. pteronyssinus
d02	D. farinae
e02/e05	Dog epithelia
e01	Cat epithelia
e03	Horse epithelia
m03	Aspergillus fumigatus
m02	Cladosporium herbarum

## Atopic 30-I

f02	Milk
f76	Alpha-Lactalbumin
f77	Beta-Lactoglobulin
f78	Casein
f01	Egg White
f75	Egg Yolk
f03	Codfish
f13	Peanut
f93	Cacao
f14	Soybean
f49	Apple
f31	Carrot
f25	Tomato
fx10	Flour Mix <sup>4</sup>
f83	Chicken Meat
fx29	Citrus Mix <sup>3</sup>
f09	Rice
gx7	6 Grass Mix <sup>1</sup>
g12	Rye pollen
e02/e05	Dog epithelia
e01	Cat epithelia
m01	Cladosporium herbarum
m06	Alternaria alternata
m03	Aspergillus fumigatus
d01	D. pteronyssinus
d02	D. farinae
t04	Hazel pollen
t03	Birch pollen
w06	Mugwort pollen
k202	CCD <sup>9</sup>

## Atopic 20-I

f02	Milk
f78	Casein
f76	Alpha-Lactalbumin
f77	Beta-Lactoglobulin
e204	Bovine serum albumin
f01/f75	Egg White / Egg Yolk Mix
f09	Rice
f14	Soybean
f92	Banana
f26	Pork meat
f27	Beef meat
f83	Chicken meat
fx10	Flour Mix <sup>4</sup>
f45	Baker's Yeast
d01/d02	D. pteronyssinus/D. farinae Mix
m02/m06	Cladosporium herbarum / Alternaria alternata Mix
t03/t07	Birch pollen / Oak pollen Mix
t02/t04	Alder pollen / Hazel pollen Mix
gx7	6 Grass Mix <sup>1</sup>
tlgE	Total IgE

## Atopic 30-II

f85	Celery
f14	Soybean
f13	Peanut
f12	Pea
f96	Avocado
f26	Pork Meat
f04	Wheat flour
f24	Shrimp
f03	Codfish
f84	Kiwi
f02	Milk
f01	Egg White
t03	Birch pollen
t09	Olive pollen
t07	White Oak pollen
t11	London Plane pollen
t23	Cypress pollen
g15	Cultivated Wheat pollen
g06	Timothy Grass pollen
w01	Ragweed pollen
w21	Wall Pellitory pollen
w06	Mugwort pollen
e01	Cat epithelia
e02/e05	Dog epithelia
m03	Aspergillus fumigatus
m06	Alternaria alternata
i06	Cockroach
d01	D. pteronyssinus
d02	D. farinae
k202	CCD <sup>9</sup>

### Korea I-Standard

f14	Soybean
f02	Milk
f01	Egg White
f23	Crab
f24/f440	Shrimp Mix <sup>5</sup>
f95	Peach
h01	House Dust
d01	D. pteronyssinus
d02	D. farinae
ix267	Cockroach Mix <sup>6</sup>
e02/e05	Dog epithelia
e01	Cat epithelia
m06	Alternaria alternata
t02/t03	Alder pollen / Birch pollen Mix
t07	White Oak pollen
w01	Ragweed pollen
w06	Mugwort pollen
g06	Timothy Grass pollen
g12	Rye pollen
tlgE	Total IgE

### Korea Standard 30

f14	Soybean
f02	Milk
f01	Egg White
f23	Crab
f24/f440	Shrimp Mix <sup>5</sup>
f95	Peach
h01	House Dust
d01	D. pteronyssinus
d02	D. farinae
ix267	Cockroach Mix <sup>6</sup>
e02/e05	Dog epithelia
e01	Cat epithelia
m06	Alternaria alternata
t02/t03	Alder pollen / Birch pollen Mix
t07	White Oak pollen
w01	Ragweed pollen
w06	Mugwort pollen
g06	Timothy Grass pollen
g12	Rye pollen
i01	Bee Venom
w22	Japanese Hop pollen
m03	Aspergillus fumigatus
m02	Cladosporium herbarum
f60	Mackerel
f11	Buckwheat flour
f13	Peanut
i81	Pupa
d72	Tyrophagus putrescentiae
k202	CCD <sup>9</sup>
tlgE	Total IgE

### Korea II-Inhalation

w07	Oxeye Daisy pollen
g09	Redtop pollen
t61	Sycamore pollen
w11	Russian Thistle pollen
t17	Japanese Cedar pollen
t215	Lilac pollen
tx251	Ash Mix <sup>7</sup>
t19	Acacia pollen
t04	Hazel pollen
w14	Pigweed pollen
g04	Meadow Fescue pollen
g02	Bermuda Grass pollen
g03	Cocksfoot pollen
w12	Goldenrod pollen
w08	Dandelion pollen
k82	Latex
m03	Aspergillus fumigatus
m02	Cladosporium herbarum
m01	Penicillium notatum
d72	Tyrophagus putrescentiae

### Korea Inhalation 30

w07	Oxeye Daisy pollen
g09	Redtop pollen
t61	Sycamore pollen
w11	Russian Thistle pollen
t215	Lilac pollen
tx251	Ash Mix <sup>7</sup>
t19	Acacia pollen
w14	Pigweed pollen
g04	Meadow Fescue pollen
g02	Bermuda Grass pollen
g03	Cocksfoot pollen
w12	Goldenrod pollen
w08	Dandelion pollen
k82	Latex
t17	Japanese Cedar pollen
d61	Tetranychus urticae
t04	Hazel pollen
t16	Pine pollen
m09	Fusarium moniliforme
d70	Acarus siro
m01	Penicillium notatum
t05	Beech pollen
t14	Cottonwood pollen
t08	Elm pollen
m05	Candida albicans
t12	Willow pollen
g08	Kentucky Blue Grass pollen
e81	Sheep epithelia
g07	Reed pollen
g01	Sweet Vernal Grass pollen

### Korea III-Food

f11	Buckwheat flour
f13	Peanut
f93	Cacao (Chocolate)
f81	Cheese, Cheddar
f60	Mackerel
f25	Tomato
f35	Potato
f40	Tuna, Yellowfin
f03	Codfish
f41	Salmon, Atlantic
f26	Pork meat
f83	Chicken meat
f27	Beef meat
fx29	Citrus Mix <sup>3</sup>
f45	Baker's Yeast
f04	Wheat flour
f09	Rice
f06	Barley flour
f47	Garlic
f299	Sweet Chestnut

### Korea Food 30

f93	Cacao (Chocolate)
f81	Cheddar
f25	Tomato
f40	Tuna, Yellowfin
f03	Codfish
f41	Salmon, Atlantic
f26	Pork meat
f83	Chicken meat
f27	Beef meat
fx29	Citrus Mix <sup>3</sup>
f45	Baker's Yeast
f04	Wheat flour
f09	Rice
f06	Barley flour
f47	Garlic
f299	Sweet Chestnut
f48	Onion
f346	Abalone
f259	Grapes
f207	Clam
f49	Apple
f281	Curry
f216	Cabbage
f313	Anchovy
f35	Potato
f244	Cucumber
f92/f91	Banana / Mango Mix
f84	Kiwi
f37	Mussel
f290	Oyster

## Mediterran 20-I

d01	D. pteronyssinus
d02	D. farinae
d71	Lepidoglyphus destructor
d72	Tyrophagus putrescentiae
e01	Cat epithelia
e02/e05	Dog epithelia
e04	Cow epithelia
e85	Chicken feathers
g12	Rye pollen
g13	Velvet Grass pollen
g21	Quack Grass pollen
g202	Maize pollen
gx8	Grass Mix*
w01	Ragweed pollen
w06	Mugwort pollen
w21	Wall Pellitory pollen
w204	Sunflower pollen
t14	Cottonwood pollen
t19	Acacia longifolia pollen
t23	Cypress pollen

## Mediterran 20-II

t02	Alder pollen
t03	Birch pollen
t04	Hazel pollen
t05	Beech pollen
t07	White Oak pollen
t09	Olive pollen
t12	Willow pollen
t16	Pine pollen
m01	Penicillium notatum
m02	Cladosporium herbarum
m03	Aspergillus fumigatus
m06	Alternaria alternata
m11	Rhizopus nigricans
gx8	Grass Mix*
h01	House Dust
i01	Bee venom
i03	Wasp venom
e81	Sheep epithelia
w06	Mugwort pollen
w07	Oxeye Daisy pollen

## Insects / CCD 6

i01	Bee venom
i03	Wasp venom
i75	Hornet venom
i71	Mosquito
i80	Black Fly
k202	CCD <sup>9</sup>

## Milk plus Gluten 6

f02	Milk
f76	Alpha-Lactalbumin
f77	Beta-Lactoglobulin
f78	Casein
e204	Bovine Serum Albumin
f79	Gluten

## rBet v1/v2 plus 20

t03	Birch pollen
t215	rBet v1
t216	rBet v2
t02	Alder pollen
t04	Hazel pollen
g06	Timothy Grass pollen
w06	Mugwort pollen
d01	D. pteronyssinus
k82	Latex
f17	Hazelnut
f13	Peanut
f85	Celery
f31	Carrot
f25	Tomato
f95	Peach
f49	Apple
f14	Soybean
f04	Wheat flour
f46	Red Pepper
fx29	Citrus Mix <sup>3</sup>

## Drugs 10

c01	Penicillin G
c02	Penicillin V
c05	Ampicillin
c06	Amoxicillin
c223	Sulfamethoxazol
c201	Cephalosporin
c436	Ofloxacin
c07	Cefactor
c103	Tetracyclin
c212	Erythromycin

## RecPollen-IT

g06	Timothy Grass pollen
g205	Phl p1
g215	Phl p5
t03	Birch pollen
t215	rBet v1
t216	rBet v2

## Explanation

<sup>1</sup> 6 Grass Mix

g03 / g04 / g05 / g06 / g08 / g13 – Cocksfoot pollen, Meadow Fescue pollen, Rye Grass pollen, Timothy Grass pollen, Kentucky Blue Grass pollen, Velvet Grass pollen

<sup>2</sup> Feathers Mix

e70 / e85 / e86 / e89 – Feathers of Goose, Chicken, Duck and Turkey

<sup>3</sup> Citrus Mix

f208 / f306 / f33 / f302 – Lemon, Lime, Orange and Mandarin

<sup>4</sup> Flour Mix

f04 / f05 / f06 / f07 – Wheat, Rye, Barley and Oat (flours that contain Gluten)

<sup>5</sup> Shrimp Mix

f24 / f440 – Atlantic and Pacific Shrimps

<sup>6</sup> Cockroach Mix

i06 / i206 / i207 – Blatella germanica, Periplaneta americana, Blatta orientalis

<sup>7</sup> Ash Mix

t25 / t15 / tx – European, White, Oregon and Black Ash

<sup>8</sup> Grass Mix

g02 / g05 / g06 – Bermuda Grass pollen, Rye Grass pollen, Timothy Grass pollen

<sup>9</sup> CCD

Cross-reactive Carbohydrate Determinants


Additional regional panels are available, please ask for further information.

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# Polycheck® Allergy Diagnostic User Manual

Scanner Based Assay for Screening and Quantification of Allergen-specific IgE

English

Allergy Panels	Panel-ID	Colour	12-Kit No.	24-Kit No.	A-Kit No.
Inhalation 10-I	11		-	05012102	05012112
Inhalation 10-II	12		-	05012202	05012212
Inhalation 10-III	48		-	05012105	05012115
Inhalation 20-I	07		-	05012002	05012022
Inhalation 30-I	36		-	05012092	05012093
Food 10-III	46		-	05012103	05012106
Food 10-IV	47		-	05012104	05012107
Food 20-I	08		-	05012001	05012021
Food 30-I	37		-	05012096	05012097
Paediatric 20-I	09		-	05012003	05012023
Atopic 20-I	31		-	05012019	05012024
Atopic 30-I	49		-	05012085	05012095
Atopic 30-II	35		-	05012090	05012091
Korea I-Standard	15		-	05012005	05012015
Korea Standard 30	54		-	05012070	05012071
Korea II-Inhalation	16		-	05012006	05012016
Korea Inhalation 30	55		-	05012072	05012073
Korea III-Food	17		-	05012007	05012017
Korea Food 30	56		-	05012074	05012075
Mediterran 20-I	33		-	05012033	-
Mediterran 20-II	34		-	05012034	-
Insects / CCD 6	50		05012314	-	05012308
Milk plus Gluten 6	52		05012318	-	05012307
rBet v1/v2 plus 20	25		-	05012009	-
Drugs 10	30		05012316	05012203	05012317
RecPollen-IT	70		05012319	-	05012320

**Packaging Unit:** 12-Kit: 1 microtiter plate with 12 tests, 24-Kit: 2 microtiter plates with 12 tests each, A-Kit: for testing with the automat

Polycheck®-Kit Components	Content	Preparation	Store at	Shelf Life
<b>Polycheck® Allergy Cassettes</b>	24 (12) Cassettes	ready to use	2 – 8 °C with desiccant in a sealed plastic bag	see expiry date
● <b>Start Solution</b> Buffered protein solution	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C	see expiry date
● <b>Anti-IgE Antibody</b> Monoclonal (murine) Antibody labelled with ligand	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C	see expiry date
○ <b>Enzyme-Labelled Anti-Ligand</b> Ligand conjugated to alkaline phosphatase	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C	see expiry date
● <b>Substrate Solution</b> 5' bromo-4' chloro-3' indolylphosphate and 4' nitroblue tetrazolium, buffered	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C protect from light	see expiry date
<b>Wash Buffer</b> Phosphate Buffer, pH 7.4	2 (1) pouches	dissolve in 1 liter demineralised water	2 – 8 °C avoid foaming	30 days after dissolving; until the expiry date, resp.

Reagents in use must be stored at 2 – 8 °C and protected from contamination. Material Safety Data Sheets are available upon request.

## Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Tests should be performed only by staff, being specially informed and trained in methods dealing with *in vitro* diagnostics. Please adhere strictly to the sequence of pipetting steps provided in this protocol. Sample material of patients (for example serum or plasma) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labelled and if necessary be handled in safety work benches (for example laminar flow bench).

## Assay Characteristics

- **Sample material:** Serum
- **Duration of assay:** Manual procedure: 2.5 hours, fully-automated procedure: about 4 hours
- **Lower limit of detection:** 0.15 kU/l
- **Specificity:** Human IgE; no cross reaction with IgG, IgA, IgD, IgM
- **Calibration:** Using single donors and pooled sera with defined allergen-specific IgE concentrations
- **Method Comparison:** Skin test and alternate solid phase immunoassays

## Materials Required

- **IT Equipment:** Personal Computer (Windows 7 or higher versions, Net Framework 3.5), printer, flatbed scanner (ask for recommended models or check on [www.polycheck.de](http://www.polycheck.de))
- **Biocheck Imaging Software (BIS):** for patient-oriented analysis, calculation and report (download from [www.polycheck.de](http://www.polycheck.de))
- **Lab Equipment:** Rocking shaker (30 rpm), adjustable pipettes for 200 up to 1000 µl; Multipette (1 ml) to dispense the wash solution
- **Deminerzalized water:** one-litre bottle for preparation of wash solution

## Specimen Preparation and Storage

Serum samples handled with professional care can be used for up to 4 weeks when stored at 2–8 °C. For longer periods samples should be stored frozen at -20 °C. Lipaemic sera should be cleared by centrifugation prior to assay.

## Method and Test Principle

The Polycheck® Screening Assay (BIS) is an enzyme immunoassay for the quantitative measurement of aller-

gen-specific IgE in serum. The relevant allergens are coated separately in lines as well as the calibrators on a carrier which is fixed in the well of the Polycheck® cassette. During incubation of the patient's serum, allergen-specific IgE binds to the corresponding allergens. Non-bound serum components are removed by washing. Monoclonal ligand-labelled anti-IgE-antibodies bind to allergen-bound IgE. Unbound antibodies are removed by washing. Enzyme-labelled anti-ligands bind to the immune complexes; surplus of enzyme conjugate is removed by washing. The substrate solution is added and specifically bound enzymes convert the colourless substrate to a dark precipitate. The colour intensity of the lines is proportional to the respective allergen-specific IgE concentration in the patient's serum. With the help of Biocheck Imaging Software (BIS) and a computer, the Polycheck® cassettes are interpreted. Each single allergen will be identified and, according to the calibrator curve present in each cassette, the concentration of each allergen-specific IgE is quantified.

## Test Performance

- ! To start the assay: all test components should be at room temperature and be mixed well.
- ! Use only reagent's lots provided with the actual kit.
- ! Powdered wash buffer has to be diluted with demineralized water at least 30 minutes prior to use. Avoid foaming.
- ! Do not allow the membranes of the test cassettes to dry during the assay.
- ! All incubation steps are performed at room temperature (18–24 °C) and with constant shaking.
- ! For the correct interpretation of the test results, a flatbed scanner with CCD sensor and a scanning resolution of 600 DPI have to be used.
- ! For test performance by automated systems additional information sheets are available.

1. Prepare a sufficient number of Polycheck® allergy cassettes and mark them – only on the long side of the cassette.
2. Moisturise the cassettes with 1 ml wash buffer, remove the wash buffer by tapping upside-down on absorbent paper.
3. Overlay allergy cassettes with 250 µl of Polycheck® start solution (blue cap) and incubate for 60 seconds (always pipet into the gap). Tap the cassettes carefully upside-down on absorbent paper.
4. Add 200 µl of the respective patient's serum into the cassette and incubate for 60 minutes on a shaker. Place the MTP-holder on the middle of the shaker

(see picture). The orientation of the MTP-holder on the shaker should remain the same during the whole test procedure.

5. Decant and wash three times with 1 ml of Polycheck® wash buffer. Tap the cassettes carefully upside-down on absorbent paper. Add 250 µl wash buffer and incubate for 5 minutes on a shaker.
6. Repeat step 5. Decant and tap the cassettes carefully.
7. Pipet 250 µl of Polycheck® anti-IgE antibody (green cap) and incubate for 45 minutes on a shaker. Decant and wash three times with 1 ml wash buffer. Tap the cassettes carefully on absorbent paper.
8. Add 250 µl Polycheck® enzyme-labelled anti-ligand (white cap) and incubate for 20 minutes on a shaker. Decant and wash as described in 7. Tap the cassettes carefully on absorbent paper.
9. Pipet 250 µl Polycheck® substrate solution (black cap) and incubate for 20 minutes in the dark. Decant and wash as described in 7.
10. Air-dry the membrane and evaluate the Polycheck® allergy cassettes using a scanner and the Biocheck Imaging Software.


for each allergen. A patient's report is generated for each patient's sample listing the allergens tested and the present concentration of allergen-specific IgE. Results are reported quantitatively in International Units (kU IgE/l). The class scores of the respective IgE concentrations and their related clinical interpretation are listed below.

**Non-specific IgE (Total IgE):** Unspecified serum-IgE concentrations may differ in a wide range in »healthy« patients. In general, for children (up to ten years) IgE levels are suspicious for atopy when exceeding 20 kU/l and in adults higher than 100 kU/l.

## Limitation of the Method

1. For a final clinical diagnosis, all clinical results and laboratory findings have to be taken into consideration.
2. Allergenicity of food allergens may be affected by cooking and/or ingestion. Consequently, IgE detection in *in vitro* diagnostic assays may be diminished or even prevented despite convincing clinical history of the patient.
3. Due to occurrence of common antigenic determinants in certain allergens of a related natural material family, the possibility of cross reactions cannot be excluded completely.

### Position of MTP-holder on Shaker

## Interpretation of Results

**Specific IgE:** For the evaluation of the Polycheck® allergy cassettes, a computer and a connected flatbed scanner have to be used. The evaluation software Biocheck Imaging Software (BIS) analyses and documents all patients' results. The MTP-holder with the Polycheck® allergy cassettes is placed upside-down on the flatbed scanner and is analysed by the software. The evaluation program checks for validity of the calibration curve, assigns the corresponding allergen to each position of the scanned band pattern and calculates the allergen-specific IgE concentration

# Interpretation of Results and Result Sheet

English

IgE [kU/l]	Class	Interpretation
< 0.35	0	No specific antibody detectable
0.35 - < 0.7	1	Very weak antibody concentration
0.7 - < 3.5	2	Weak antibody concentration
3.5 - < 17.5	3	Clear antibody concentration
17.5 - < 50	4	Strong antibody concentration
50 - < 100	5	Very strong antibody concentration
≥ 100	6	Extremely high antibody concentration

Dr. Jana Smith  
Laboratory  
Any Street  
12345 Any City

Quantitative Multiparameter Assay  
ISO 13485:2012 certified  
www.polycheck.de

Name: Doe  
First Name: John  
ID or d.o.b.: 15.10.1972  
File Name: ...nbankn/Bespiel Inha Flyer.bio  
Scan of 02.11.2015, class. # 1 (V4.01.02)


Test (Allergy): Inhalation 20-I  
Order # / Origin: 12345 /  
Date of Order: 02.11.2015  
Print Date: 26.09.2016 (V. 4.04.12)  
M:\Allergien\neue Panels\Korea-Panel\OCV\Bocheck\Btmapa\02\_11\_2015  
113\_48\_16\_01.bmp

Allergens	Class	IgE [kU/l]	Class				
			0	1	2	3	4+5
t03 Birch pollen	5	91	[Bar chart showing Class 5]				
t02 Alder pollen	5	63	[Bar chart showing Class 5]				
t04 Hazel pollen	3	11	[Bar chart showing Class 3]				
t07 Oak pollen	0	0.16	[Bar chart showing Class 0]				
g06 Timothy Grass pollen	1	0.48	[Bar chart showing Class 1]				
g12 Rye pollen	1	0.49	[Bar chart showing Class 1]				
w06 Mugwort pollen	2	0.88	[Bar chart showing Class 2]				
w09 Plantain pollen	0	<0.15	[Bar chart showing Class 0]				
d01 D. pteronyssinus	0	<0.15	[Bar chart showing Class 0]				
d02 D. farinae	0	<0.15	[Bar chart showing Class 0]				
e02/e05 Dog epithelia	0	<0.15	[Bar chart showing Class 0]				
e01 Cat epithelia	0	<0.15	[Bar chart showing Class 0]				
e03 Horse epithelia	4	24	[Bar chart showing Class 4]				
e06 Guinea pig epithelia	1	0.50	[Bar chart showing Class 1]				
e84 Hamster epithelia	0	<0.15	[Bar chart showing Class 0]				
e82 Rabbit epithelia	0	<0.15	[Bar chart showing Class 0]				
m03 Asp. fumigatus	2	1.2	[Bar chart showing Class 2]				
m02 Cladosp. herbarum	0	<0.15	[Bar chart showing Class 0]				
m01 Pen. notatum	0	<0.15	[Bar chart showing Class 0]				
m06 Alt. alternata	0	<0.15	[Bar chart showing Class 0]				

Class	Conc. IgE [kU/l]	Explanation	Class	Conc. IgE [kU/l]	Explanation
0	<0.35	No specific antibody detection	4	17.5 - <50	Strong antibody concentration
1	0.35 - <0.7	Very weak antibody concentration	5	50 - <100	Very strong antibody concentration
2	0.7 - <3.5	Weak antibody concentration	6	≥ 100	Extremely high antibody concentration
3	3.5 - <17.5	Clear antibody concentration			

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

Allergene	Klasse	IgE [kU/l]	Klasse				
			0	1	2	3	4 + 5
t03 Birke	5	91	[Bar chart showing high concentration]				
t02 Erle	5	63	[Bar chart showing high concentration]				
t04 Hasel	3	11	[Bar chart showing moderate concentration]				
t07 Eiche	0	0.16	[Bar chart showing very low concentration]				
g06 Lieschgras	1	0.48	[Bar chart showing very low concentration]				
g12 Roggen	1	0.49	[Bar chart showing very low concentration]				
w06 Beifuss	2	0.88	[Bar chart showing low concentration]				
w09 Wegerich	0	<0.15	[Bar chart showing very low concentration]				
d01 D.pteronissinus	0	<0.15	[Bar chart showing very low concentration]				
d02 D.farinae	0	<0.15	[Bar chart showing very low concentration]				
e02/e05 Hundeseithelien	0	<0.15	[Bar chart showing very low concentration]				
e01 Katzenepithelien	0	<0.15	[Bar chart showing very low concentration]				
e03 Pferdeepithelien	4	24	[Bar chart showing high concentration]				
e06 Meerschweinchenepith.	1	0.50	[Bar chart showing very low concentration]				
e84 Hamsterepithelien	0	<0.15	[Bar chart showing very low concentration]				
e82 Kaninchenepithelien	0	<0.15	[Bar chart showing very low concentration]				
m03 Asp.fumigatus	2	1.2	[Bar chart showing low concentration]				
m02 Cladosp.herbarum	0	<0.15	[Bar chart showing very low concentration]				
m01 Pen.notatum	0	<0.15	[Bar chart showing very low concentration]				
m06 Alt.alternata	0	<0.15	[Bar chart showing very low concentration]				

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Scanne de 02.11.2015, cass. # 1 (V4.01.02)

Allergènes	Classe	IgE [kU/l]	Classe				
			0	1	2	3	4+5
t03 Bouleau	5	91	[Bar chart showing high concentration]				
t02 Auline	5	63	[Bar chart showing high concentration]				
t04 Noisetier	3	11	[Bar chart showing low concentration]				
t07 Chêne	0	0.16	[Bar chart showing no concentration]				
g06 Phléole	1	0.48	[Bar chart showing very low concentration]				
g12 Seigle	1	0.49	[Bar chart showing very low concentration]				
w06 Armoise	2	0.88	[Bar chart showing low concentration]				
w09 Plantain	0	<0.15	[Bar chart showing no concentration]				
d01 D.pteronysinus	0	<0.15	[Bar chart showing no concentration]				
d02 D.farinae	0	<0.15	[Bar chart showing no concentration]				
e02/e05 Chien	0	<0.15	[Bar chart showing no concentration]				
e01 Chat	0	<0.15	[Bar chart showing no concentration]				
e03 Cheval	4	24	[Bar chart showing high concentration]				
e06 Cobaye	1	0.50	[Bar chart showing very low concentration]				
e84 Hamster	0	<0.15	[Bar chart showing no concentration]				
e82 Lapin	0	<0.15	[Bar chart showing no concentration]				
m03 Asp.fumigatus	2	1.2	[Bar chart showing low concentration]				
m02 Cladosp.herbarum	0	<0.15	[Bar chart showing no concentration]				
m01 Pen.notatum	0	<0.15	[Bar chart showing no concentration]				
m06 Alt.alternata	0	<0.15	[Bar chart showing no concentration]				

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**Polycheck® is a Brand of:** Biocheck GmbH / Vorbergweg 41 / 48159 Münster / Germany / **Version: 010617**