Polycheck® Autoimmune Diagnostics Multiparameter Technology

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Order Numbers

Our Packaging Units

12-Kit: 1 microtiter plate with 12 tests each
24-Kit: 2 microtiter plates with 12 tests each
A-Kit: for testing with the automat
Semi-A.-Kit: for testing with the semi-automat

Autoimmune Kits	Colour	ID	12-Kit No.	24-Kit No.	A-Kit No.	Semi-A. No.
AI-STD		40	05012050	05012040	05012060	05012560
AI-STD-X		44	05012054	05012044	05012064	05012564
Rheuma		53	05012056	05012046	05012066	05012566
Liver/Kidney		42	05012052	05012042	05012062	05012562
Myositis		41	05012051	05012041	05012061	05012561
Celiac IgG		43	05012053	05012043	05012063	05012563
Celiac IgA +total		74	05012057	05012047	05012067	05012567
Celiac Dual Kit*		74+43	-	05012247	05012267	05012568

*Celiac Dual Kit: 12 tests Celiac IgA +total and 12 tests Celiac IgG

Autoimmune Panels

AI-STD

01	Ro/SS-A 60
02	La/SS-B
03	Scl-70
04	PM/Scl 100
05	Sm
06	PCNA
07	dsDNA
08	Rib. Phosphopr. Po
09	CENP-B
10	AMA M2
11	PR3
12	MPO
13	TPO
14	Thyroglobulin [TG]

Rheuma

01	Ro/SS-A 52
02	Jo-1
03	Sm
04	CENP-B
05	La/SS-B
06	Histone
07	PM/Scl 100
08	Rib. Phosphopr. Po
09	dsDNA
10	Scl-70
11	U1-snRNP
12	Ro/SS-A 60
13	Mi-2
14	Ku

AI-STD-X

01	Ro/SS-A 52
02	La/SS-B
03	CENP-B
04	Scl-70
05	dsDNA
06	Jo-1
07	MPO
08	PR3
09	AMA M2
10	LC 1
11	LKM 1
12	PM/Scl 100
13	SRP 54
14	SP 100
15	gp 210
16	Ku
17	Sm
18	U1-snRNP

Liver/Kidney

01	AMA M2
02	SP 100
03	gp 210
04	LKM 1
05	LC 1
06	CENP-B
07	GBM
08	PR3
09	MPO
10	SLA/LP

Myositis

01	Mi-2
02	PM/Scl 100
03	Jo-1
04	SRP 54
05	PL-7
06	PL-12
07	Ku

Celiac IgG

01	deam. Gliadin IgG
02	tTG IgG
03	Intrinsic Factor

Celiac IgA +total

01	deam. Gliadin IgA
02	tTG IgA
03	total IgA

Celiac Dual Kit

12 x	Celiac IgA +total
12 x	Celiac IgG

Polycheck® Autoimmune Diagnostic User Manual

Scanner Based Assay for Screening and Quantification of Antigen-specific Antibodies

Polycheck®-Kit Components	Content	Preparation	Store at	Shelf Life
Polycheck® Autoimmune Cassettes	24 [12] Cassettes	ready to use	2 – 8 °C with desiccant in a sealed plastic bag	see expiry date
Start Solution Buffered protein solution	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C	see expiry date
Anti-IgG-Antibody Anti-IgA-Antibody Alkaline phosphatase conjugated detection antibody	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C	see expiry date
Substrate Solution 5'bromo-4'chloro-3' indolylphosphate and 4' nitroblue tetrazolium, buffered	Manual: 2 (1) x 3.5 ml Automat: 1 x 10 ml	ready to use	2 – 8 °C protect from light	see expiry date
Diluent Buffered protein solution	Manual: 2 (1) x 20ml Automat: 1 x 20 ml	ready to use	2 – 8 °C	see expiry date
Wash Buffer Phosphate Buffer, pH 7.4	Manual: 2 (1) pouches Automat: 1 pouch	dissolve in 1 litre demineralised water	2 – 8 °C avoid foaming	30 days after dissolving; until the expiry date, resp.

Reagents in use must be stored at 2 – 8 °C and protected from contamination. Material Safety Data Sheets are available upon request.

Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Tests should be performed only by staff, being specially informed and trained in methods dealing with *in vitro* diagnostics. Please adhere strictly to the sequence of pipetting steps provided in this protocol. Sample material of patients (for example serum or plasma) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labelled and if necessary be handled in safety work benches (e.g. laminar flow bench).

Assay Characteristics

- **Sample material:** Serum or plasma (EDTA, heparin, citrate), dilution 1:100
- **Duration of assay:**
Manual procedure: 1 hour and 45 minutes
Fully-automated procedure: about 3.5 hours
- **Lower limit of detection:** 0.15 kU/l
- **Specificity:** Specific for human IgG; respectively specific for human IgA (in case of Celiac IgA +total panel)
- **Calibration:** Using single donors and pooled sera with defined antigen-specific IgG or IgA concentrations, respectively
- **Method Comparison:** indirect immunofluorescence tests and alternate solid phase immunoassays

Materials Required

- **IT Equipment:** Personal Computer (Windows 7 or higher versions, Net Framework 3.5), printer, flatbed scanner (ask for recommended models or check on www.polycheck.de)
- **Biocheck Imaging Software (BIS):** for patient-oriented analysis, calculation and report (download from www.polycheck.de)
- **Lab Equipment:** Rocking shaker (30 rpm), adjustable pipettes for 5 up to 1000 µl; Multipipette (1 ml) to dispense the wash solution
- **Deminerlized water:** one-litre bottle for preparation of wash solution

Specimen Preparation and Storage

Serum samples handled with professional care can be used for up to 4 weeks when stored at 2–8 °C. For longer periods samples should be stored frozen at -20 °C. Lipaemic sera should be cleared by centrifugation prior to assay. For autoimmune diagnostics all sera have to be diluted 1:100 with dilution buffer before using them in the assay.

Method and Test Principle

The Polycheck® Autoimmune Screening Assay is an enzyme immunoassay for the quantitative measurement of antigen-specific IgG or IgA in serum. The relevant antigens as well as the calibrators are coated separately in lines on a carrier which is fixed in the well of the Polycheck® cassette. During incubation of the patient's serum, antigen-specific IgG or IgA binds to the corresponding antigens. Non-bound serum components are removed by washing. Subsequently, alkaline phosphatase conjugated detection antibodies identify antigen-bound IgG or IgA, respectively. Unbound antibodies are removed by washing. The substrate solution is added and specifically bound enzymes convert the colourless substrate to a dark precipitate. The colour intensity of the lines is proportional to the respective antigen-specific IgG or IgA concentration in the patient's serum. With the help of Biocheck Imaging Software (BIS), a scanner and a computer, the Polycheck® cassettes are interpreted. Each single antigen will be identified and, according to the calibrator curve present in each cassette, the concentration of each antigen-specific IgG or IgA is quantified.

Test Performance


- ! To start the assay: all test components should be at room temperature and be mixed well.
- ! Use only reagent's lots provided with the actual kit.
- ! Powdered wash buffer has to be diluted with demineralized water at least 30 minutes prior to use. Avoid foaming.
- ! Do not allow the membranes of the test cassettes to dry during the assay.
- ! All incubation steps are performed at room temperature (18–24 °C) and with constant shaking.
- ! For the correct interpretation of the results, a flatbed scanner with CCD sensor and a scanning resolution of 600 DPI have to be used.
- ! For test performance by automated systems additional information sheets are available.

1. Prepare a sufficient number of Polycheck® autoimmune cassettes and mark them – only on the long side of the cassette.
2. Moisturise the cassettes with 1 ml wash buffer, remove the wash buffer by tapping upside-down on absorbent paper.
3. Overlay autoimmune cassettes with 250 µl of Polycheck® start solution (blue cap) and incubate for 60 seconds (always pipet into the gap, see

scheme). Tap the cassettes carefully upside-down on absorbent paper.

4. Add 250 µl of the respective diluted patient's serum – **dilution: 1:100, e.g. 15 µl serum + 1485 µl Polycheck® diluent (lilac cap)** – into the cassette and incubate for 45 minutes on a shaker. Place the MTP-holder on the middle of the shaker (see picture). The orientation of the MTP-holder on the shaker should remain the same during the whole test procedure.
5. Decant and wash three times with 1 ml of Polycheck® wash buffer. Tap the cassettes carefully upside-down on absorbent paper. Add 250 µl wash buffer and incubate for 5 minutes on a shaker.
6. Repeat step 5. Decant and tap the cassettes carefully.
7. Wash three times with 1 ml of Polycheck® wash buffer. Tap the cassettes carefully upside-down on absorbent paper.
8. Pipet 250 µl of Polycheck® detection antibody (Anti-IgG-Antibody: red cap or Anti-IgA-Antibody: yellow cap) and incubate for 30 minutes on a shaker. Decant and wash three times with 1 ml wash buffer. Tap the cassettes carefully on absorbent paper.
9. Pipet 250 µl Polycheck® substrate solution (black cap) and incubate for 20 minutes in the dark. Decant and wash as described in step 7. Tap the cassettes carefully on absorbent paper.
10. Airdry the membrane for at least 30 min and evaluate the Polycheck® autoimmune cassettes using a scanner and the Biocheck Imaging Software.

Position of MTP-holder on Shaker

Interpretation of Results

For the evaluation of the Polycheck® autoimmune cassettes, a computer and a connected flatbed scanner have to be used. The evaluation software Biocheck Imaging Software (BIS) analyses and documents all patients' results. The MTP-holder with the Polycheck® autoimmune cassettes is placed upside-down on the flatbed scanner and is analysed by the software. The evaluation program checks for validity of the calibration curve, assigns the corresponding antigen to each position of the scanned band pattern and calculates the antigen-specific IgG or IgA concentration for each antigen. A patient's report is generated for each patient's sample listing the antigens tested and the present concentration of antigen-specific IgG or IgA, respectively. Results are reported quantitatively in International Units (kU IgG/l or kU IgA/l, respectively). The range scores of the respective IgG or IgA concentrations and their related interpretation are listed below.

Additional Information

1. For a final clinical diagnosis, all clinical results and laboratory findings have to be taken into consideration.
2. Some antibodies appear years before outbreak of the respective disease. Therefore, detection of these antibodies can be a significant indicator for the onset of a disease.

For *total IgA*, other definitions of ranges apply. Values > 0.8 kU/l are interpreted as normal IgA response.

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Polycheck® is a Brand of: Biocheck GmbH / ☎ +49 251 21 50 869 / ✉ office@polycheck.de
Polycheck® is a Brand of Biocheck GmbH
Biocheck GmbH is certified to EN ISO 13485:2016